

Accessories



Feedthrough

Feedthrough assembly to pass 20/30 thermocouples into autoclaves

Order Reference: **PRSET20/ PRSET30**



Shelf for Loggers Calibration

Shelf for data loggers calibration in dry blocks

Order Reference: **DLS**



Y/DN50

Stainless steel Y adapter. Connects:
- 1 port to the autoclave
- 1 port to a feedthrough
- 1 port to the pressure transducer

Order Reference: **Y/DN50**



Pressure Transducer

EN285 compliant triclamp pressure transducer to connect on the autoclave during validations

Order Reference: **PBMH25**



Dry Block Insert

Specific insert for validation in dry blocks. Insert 20 Loggers/Tc & one reference probe

Order Reference: **AL, BR**



Shelf Mapping Cap

Temperature transfer adapter from the shelf to the logger's sensor during freeze dryers mappings

Order Reference: **LyoPuck_RF, LyoPuck_ST**



Fedegari Adapter

Adapter for Fedegari autoclave to feedthrough

Order Reference: **AD-PREST20**



Pliers

Tools to open loggers for battery replacement

Order Reference: **PL-ST, PL-RF, PL-Bio, PL-LMS**



Bracket

Mounting bracket for data loggers

Order Reference: **BR-RF, BR-ST**



Repeater

Repeater to extend the coverage area of the Access Point

Order Reference: **REX**



GMP Pipe

Adapter to insert data logger's sensor inside pipes

Order Reference: **GMP_1.5_PIPE**
GMP_2_PIPE



XpertVal Module

Module for connecting thermocouples to XpertVal data logger

Order Reference: **2638A-100**



Activation Point

Activates GMP RF and LMS wireless dataloggers

Order Reference: **ActP**



Lyo Caterpillar

Allows to measure product temperature inside a freeze-dryer (lyophilizer) and eliminate the influence from data logger

Order Reference: **VCX-XXR-XXml**



Sensor Holder

Cap with holder for flexible temperature sensor inside the vial
xx = vial neck diameter (mm)
yyy = sensor support length

Order Reference: **CSS-xx-yyy**



Twin Borell

Cap holders for product measurement inside vials

Order Reference: **TW BC_XXXml**